Exploring the Effects of Indium in a Four Element Magnetic Alloy Kenny Au, John Martin, Dr. Peter Mullner

Why? To make magnets CHEER. Cheaper, Helpful, Efficient, **Environment-Friendly, Reliable**



This study addresses the growing concern with the difficulty of obtaining Nd for NdFeB magnets. Mn-Al has the potential to be a substitute but it lacks the high coercivity and ferromagnetic phase stability. Indium and gallium help with coercivity and phase stability, respectively. Through a two-step heat annealing process, the quaternary alloy showed an increase in coercivity although a greater decrease in magnetization when compared to a sample of the same composition without indium.

Background





- NdFeB demand will greatly increase; non rare-earth elements can be used to create novel magnets
- The magnetization depends on Mn atom magnetic moments being aligned - found only in this crystal structure
- · Ga and In can substitute for AI because of their



0.25

0.2

0.1

n

N



EDS Linescan across boundary precipitate

Liu

MCR

AR 277

Coercivity Plot

- The best order of heat treatments along with cooling methods is 1100°C quenched then 450 °C air
- cooled. Two heat treatments need to be done to achieve the highest coercivity
- Time plays a role as a 4 hour heat treatment starts to show ferromagnetism while a 24 hour
- treatment decreases the coercivity.



Scanning Electron Microscope

BOISE STATE UNIVERSITY



Through EDS analysis, the circular precipitate has a very high concentration of indium. This suggests that these precipitates act as a pin in the alloy thus increasing the coercivity of the alloy.

Magnetization Comparisons



- Conclusion:
- Indium increases the coercivity of the alloy 0.2 T to 0.06T in a non-indium alloy sample. However the magnetization is very low.
- 2 Heat treatments need to be done in a specific order in order to maximize the coercivity and the magnetization of the alloy

Future Works

- If the composition is the same then test different temperatures with different times
- If composition changes then repeat the best heat treatment course from here and compare

Acknowledgement: This work was supported by the NSF via the

Processing and Results

Sample AR 300E, a Mn-Al-Ga-In alloy (55-36-6-3 at%) micrograph taken at 10X zoom.